




Puritan®

Quality since 1919





Puritan®
MEDICAL

Healthcare professionals have been relying on Puritan® since 1919. And our product line has grown over that time, evolving to meet the changes in the medical field.

This line includes everything from tongue depressors, still made with fine, splinter-free Northern White Birch, to our convenient wound measuring device, and other specialty products.

We will continue to develop and offer standard and custom products backed by our reputation for quality and consistency.

The broadest single-use line available

Where it all started – quality, consistency, and reliability.

Our broad medical line includes traditional products we've been manufacturing for decades, market-driven innovations such as our foam wound measuring device, and everything in between.

US MADE

We card and spin tip fibers in-house.

CONVENIENCE

Packaging designed for easy identification and inventory control.



CONSISTENCY

We guarantee reliable products.

FEATURED PRODUCTS:

Left: Cotton-tipped applicator (#806-WC packaging)









Right: Sterile cotton-tipped applicator (#25-806 1WC) includes one swab inside a sterile wrapper

COTTON-TIPPED APPLICATORS

Available in a variety of handle materials, sizes and configurations. For these applicators, we only use US-grown cotton. The cotton is carded, spun and applied in our facility to assure you of a consistent, high-quality product, order after order.





	Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	non-sterile sterile				
	800 C	Mini tip, 0.035" diameter aluminum handle	17.45 x 1.981 mm 0.687 x 0.078"	165.1 mm 6.547"	10/100, 6 lbs
	25-800 C 50	Mini tip, 0.035" diameter aluminum handle, curled end	17.45 x 1.981 mm 0.687 x 0.078"	138.1 mm 5.437"	10/50/1, 4 lbs
	803-PC	Standard tip, polystyrene handle	17.45 x 4.369 mm 0.687 x 0.172"	79.375 mm 3.125"	10/10/50, 7 lbs
	25-803 2PC				
	803-PCL	Large tip, polystyrene handle	17.449 x 6.35 mm 0.687 x 0.25"	76.2 mm 3"	10/5/50, 5 lbs
	803-WC	Standard tip, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	76.2 mm 3"	10/10/100, 10 lbs
	25-803 2WC				10/100/2, 5 lbs
	803-WCL	Large tip, wood handle	17.449 x 5.943 mm 0.687 x 0.234"	76.2 mm 3"	10/5/100, 7 lbs
	804-PC UNBRK	Standard tip, unbreakable polystyrene handle	17.449 x 5.156 mm 0.687 x 0.203"	101.6 mm 4"	10/10/100, 12 lbs
	806-PC	Standard tip, polystyrene handle	17.449 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"	10/10/100, 17 lbs
	25-806 1PC				10/100/1, 6 lbs
	25-806 2PC				10/100/2, 7 lbs
	806-PCL	Large tip, polystyrene handle	17.449 x 6.35 mm 0.687 x 0.25"	152.4 mm 6"	10/5/100, 10 lbs
	806-WC	Standard tip, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"	10/10/100, 16 lbs
	25-806 1WC				10/100/1, 6 lbs
	25-806 2WC				10/100/2, 7 lbs
	25-806 5WC				10/20/5, 5 lbs
	25-806 10WC				10/20/10, 6 lbs

COTTON-TIPPED APPLICATORS, CONTINUED

	Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	non-sterile sterile				
	806-WCL	Large tip, wood handle	17.449 x 5.943 mm 0.687 x 0.234"	152.4 mm 6"	10/5/100, 10 lbs
	25-806 10WCL				10/20/10, 6 lbs
	25-806 WCLB 10/50				10/50, 2 lbs
	806-WCXL	Extra-large tip, wood handle	17.449 x 7.924 mm 0.687 x 0.312"	150.79 mm 5.937"	10/10/50, 12 lbs
	821-WC DBL	Double-tipped (one standard and one tapered), wood handle	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"	10/5/100, 9 lbs
	823-WC	Tapered mini tip, wood handle	12.7 x 3.048 mm 0.5 x 0.12"	76.2 mm 3"	10/10/100, 8 lbs
	826-WC	Tapered mini tip, wood handle	12.7 x 3.048 mm 0.5 x 0.12"	152.4 mm 6"	10/10/100, 12 lbs
	826-WC DBL	Double tapered mini tip, wood handle	12.7 x 3.048 mm 0.5 x 0.12"	152.4 mm 6"	10/10/100, 12 lbs
	833-WCS	Small tip, wood handle, in plastic vial	15.875 x 3.962 mm 0.625 x 0.156"	76.2 mm 3"	50/100, 8 lbs
	869-WCS	Extra-small tip, wood handle	12.7 x 3.175 mm 0.5 x 0.125"	152.4 mm 6"	10/10/100, 13 lbs

OB/GYN AND PROCTOSCOPIC APPLICATORS

Oversized and extra-absorbent, these swabs are perfect for OB/GYN, proctoscopic, veterinary, and forensic use.

	Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	non-sterile sterile				
	808	OB/GYN applicator, rayon tip, paper handle	28.575 x 12.7 mm 1.125 x 0.5"	207.975 mm 8.188"	10/50, 5 lbs
	808 100				10/100, 8 lbs
	808 BULK				500, 4 lbs
	808 COTTON	OB/GYN applicator, cotton tip, paper handle	28.575 x 12.7 mm 1.125 x 0.5"	207.975 mm 8.188"	10/50, 5 lbs
	25-808 1PR	OB/GYN applicator, rayon tip, polystyrene handle	28.575 x 12.7 mm 1.125 x 0.5"	203.2 mm 8"	10/50/1, 8 lbs
	25-808 2PR				10/50/2, 11 lbs
	816	Proctoscopic applicator, rayon tip, polystyrene handle	28.575 x 12.7 mm 1.125 x 0.5"	406.4 mm 16"	10/50, 8 lbs
	816 100				10/100, 13 lbs
	816 BULK				500, 6 lbs

CYTOLOGY COLLECTION DEVICES

Puritan® cytology products are made of the very best materials, in the most effective shapes, to maximize sample collection and provide a comfortable patient experience at the same time.


	Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	non-sterile sterile				
	2188	Histobrush®, soft nylon bristles, tapered tip	19.05 x 5.08 mm 0.75 x 0.2"	203.2 mm 8"	10/100, 5 lbs
	25-2188				10/100/1, 8 lbs
	2195	Rovers® Cervex-Brush®	28.575 x 19.05 mm 1.25 x 0.75"	203.2 mm 8"	35/25, 6 lbs
	2195 COMBI	Rovers® Cervex-Brush® ecto and endo	34.29 x 19.05 mm 1.35 x 0.75"	203.2 mm 8"	4/25, 5 lbs
	2196	Rovers® EndoCervex-Brush®	33.782 x 5.435 mm 1.333 x 0.214"	203.2 mm 8"	10/50, 3 lbs
	2197	CerviSoft®, soft foam tip	25.4 x 4.775 mm 1 x 0.188"	179.374 mm 7.062"	10/100, 3 lbs
	25-2197				10/100/1, 7 lbs
	2199	Histobrush®, soft nylon bristles, tapered tip	19.05 x 5.08 mm 0.75 x 0.2"	177.8 mm 7"	10/100, 4 lbs
	25-2199				10/100/1, 8 lbs

FOAM WOUND MEASURING DEVICE (DM STICK)

The soft foam-tipped applicator is flexible to conform to wound shape and features a metric rule for easy measuring. Designed to assess decubitus ulcers but also perfect for forensics, emergency rooms, long-term care, and more.


FEATURED PRODUCT: Foam wound measuring device (#1506-PFB DM)



	Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	non-sterile sterile				
	1506-PFB DM	DM stick, foam tip, polystyrene handle with printed graduations	9.98 x 5.08 mm 0.393 x 0.2"	152.4 mm 6"	20/50, 2 lbs
	25-1506 IPF DM				4/50/1, 2 lbs




WOODEN SCRAPERS

These scrapers may be used for specimen collection and treating lesions in areas such as the throat or skin.

	Item No.	Product Description	Full Length	Packaging pack, weight
non-sterile sterile				
	2173	Ayres style	177.8 mm 7"	10/500, 23 lbs
	25-2173			10/100/1, 9 lbs






WOODEN APPLICATOR STICKS

Made for general-purpose applications. Choose the size and shape for your specimen transfer, mixing rods, and laboratory work.

	Item No.	Product Description	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile				
	805	Wood	304.8 x 2.108 mm 12 x 0.083"	12/1000, 23 lbs
	807	Wood	147.624 x 2.108 mm 5.182 x 0.083"	20/1000, 18 lbs
	25-807 2W			10/100/2, 6 lbs
	807 HOSPITAL	Wood, straight cut ends, economy-grade, HOSPITAL brand		20/1000, 19 lbs
	809	Wood	304.8 x 3.175 mm 12 x 0.125"	12/500, 25 lbs


WOODEN CUTICLE/ORANGE STICKS

Ideal as an all-purpose stick for most clinical, diagnostic and forensic settings.
Excellent for podiatric practices and use in the beauty industry.

	Item No.	Product Description	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile				
	2910	Wood, double bevel both ends	177.8 x 3.962 mm 7 x 0.156"	10/100, 5 lbs
	2911	Wood, double bevel both ends, wrapped		4/125/1, 6 lbs
	2925	Wood, one single bevel end, one pointed end	152.4 x 3.962 mm 6 x 0.156"	10/100, 4 lbs
	2930	Wood, one single bevel end, one double bevel end	177.8 x 3.962 mm 7 x 0.156"	10/100, 5 lbs
	2950	Wood, double bevel both ends	177.8 x 4.775 mm 7 x 0.188"	10/100, 6 lbs
	2970	Wood, one single bevel end, one double bevel end	177.8 x 4.775 mm 7 x 0.188"	10/100, 6 lbs

WOODEN TRANSFER DEVICE

A simple medical device to replace the common toothpick. Non-sterile or convenient sterile multi-pack for manual colony transfer for MALDI-TOF and other applications.

	Item No.	Product Description	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile				
	28107 PTD 5/600	Pointed tip	73 x 2.108 mm 2.875 x 0.083"	5/600, 2 lbs
	25-28107			10/40/26, 2 lbs

WOODEN TONGUE DEPRESSORS

Northern White Birch produces a smooth, splinter-free finish, has no taste or odor, and maintains a high tensile strength. Used for a wide variety of applications such as patient oral exams, mixing or applying medicines, crafting and/or industrial use.

	Item No.	Product Description	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile				
	700	Infant	114.3 x 9.53 mm 4.5 x 0.375"	10/1000, 37 lbs
	25-700			10/100/1, 8 lbs
	704	Standard	152.4 x 17.475 mm 6 x 0.688"	10/500, 35 lbs
	704 BULK			3,700, 28 lbs (bulk)
	704 BRIGHTWOOD	Standard, economy-grade, Brightwood® brand	152.4 x 17.475 mm 6 x 0.688"	10/500, 35 lbs
	704 HOSPITAL	Standard, economy-grade, HOSPITAL brand	152.4 x 17.475 mm 6 x 0.688"	10/500, 35 lbs
	705	Regular	152.4 x 19.05 mm 6 x 0.75"	10/500, 38 lbs
	25-705			10/100/1, 10 lbs
	706	Standard	152.4 x 17.475 mm 6 x 0.688"	12/100, 9 lbs
	709	Standard, individually wrapped	152.4 x 17.475 mm 6 x 0.688"	10/250, 21 lbs
	710	Junior, individually wrapped	139.7 x 15.875 mm 5.5 x 0.625"	10/250, 21 lbs
	710 CHERRY	Junior, individually wrapped, flavored		20/50, 10 lbs
	710 GRAPE			
	710 STRAWBERRY			
	711	Junior	139.7 x 15.875 mm 5.5 x 0.625"	10/500, 29 lbs
	25-711			10/100/1, 10 lbs
	712	Junior	139.7 x 15.875 mm 5.5 x 0.625"	12/100, 8 lbs

Did you know?

One average-size birch tree
yields 80,000 ice cream sticks,
40,000 tongue depressors, or
200,000 6" applicator sticks.



FEATURED PRODUCTS:

Left to right: 6" Standard wooden
tongue depressor (#706); 5.5"
Junior wooden tongue depressor
(#711); Infant wooden 4.5"
tongue depressor (#700)

Puritan Medical Products Co.
31 School Street
PO Box 149
Guilford, Maine 04443-0149 USA

800-321-2313 (US & Canada)
207-876-3311
puritanmedproducts.com
sales@puritanmedproducts.com